

Blank-SAP

NON-TARGETED SAPORIN CONTROL MOLECULE

*a tool for use as control for peptide-targeted lesioning agents;
non-targeted via non-specific peptide conjugated to saporin*

Catalog Number: IT-21
Quantity: 25 micrograms, 100 micrograms, 250 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative. Sterile-filtered.

Background:

Controls are a vital part of the scientific procedure; without them it is difficult to isolate the specific effects from the non-specific or artifactual. This control molecule is the same molecular weight, consists of similar, comparable materials and is synthesized with the same protocols as the targeted conjugates. The difference is the cell-specific targeting agents are replaced with "blanks," antibodies or peptides that have no specificity, and no ability to target cells. In short, they are the perfect control molecules for behavioral experiments with Advanced Targeting Systems' targeted conjugates.

Blank-SAP is the perfect control for use with peptide conjugates. The sequence of the non-targeted peptide of this molecule is an 11-amino acid, randomly mixed version of the sequence of melanocyte-stimulating hormone, with amino acid residues that are typical of peptides that bind to G-protein-coupled receptors. Examination of the peptide sequence using NCBI/BLAST reveals no homologous sequences. Thus, Blank-SAP can be used as control for any non-specific effects of the toxin and provides a definitive baseline for determining the effects of a targeted peptide conjugate.

Specificity and Preparation:

This control conjugate (molecular weight 32 kDa) has no known specificity. Blank-SAP is a chemical conjugate between a non-targeted peptide and the ribosome-inactivating protein, saporin. The product is routinely tested by cytotoxicity assay.

Usage and Storage:

Blank-SAP serves as a control for peptide-targeted toxins (SSP-SAP, dermorphin-SAP / MOR-SAP, CRF-SAP, NPY-SAP, CCK-SAP, Bombesin-SAP, Oxytocin-SAP, and Galanin-SAP). **There may be lot-to-lot variation in material; working dilutions must be determined by end user. If this is a new lot, you must assess the proper working dilution before beginning a full experimental protocol.**

Gently spin down material before use; 5-10 seconds in a microfuge should be adequate. The material should be stored at -20°C for one year. Avoid repeated freezing and thawing.

Do not use a reducing agent (such as dithiothreitol, beta-mercaptoethanol or ascorbic acid) with this material. It will inactivate the toxin.

For disposal: autoclave, or expose to 0.2 M NaOH, materials that come into contact with the toxin.

Blank-SAP
NON-TARGETED SAPORIN CONTROL MOLECULE

References:

1. Tazzari PL, Bolognesi A, De Toter D, Falini B, Lemoli RM, Soria MR, Pileri S, Gobbi M, Stein H, Flenghi L, *et al.* (1992) Ber-H2 (anti-CD30)-saporin immunotoxin: a new tool for the treatment of Hodgkin's disease and CD30+ lymphoma: *in vitro* evaluation. *Brit J Haemat* 81:203-211.
2. Dinota A, Tazzari PL, Michieli M, Visani G, Gobbi M, Bontadini A, Tassi C, Fanin R, Damiani D, Grandi M, *et al.* (1990) *In vitro* bone marrow purging of multidrug-resistant cells with a mouse monoclonal antibody directed against Mr 170,000 glycoprotein and a saporin-conjugated anti-mouse antibody. *Cancer Res* 50:4291-4294.
3. Thorpe PE, Brown AN, Bremner JA Jr, Foxwell BM, Stirpe F (1985) An immunotoxin composed of monoclonal anti-Thy 1.1 antibody and a ribosome-inactivating protein from *Saponaria officinalis*: potent antitumor effects *in vitro* and *in vivo*. *J Natl Cancer Inst* 75(1):151-159.


Safety:

Good laboratory technique must be employed for safe handling of this product.

This requires observation of the following practices:

1. Wear appropriate laboratory attire, including lab coat, gloves and safety glasses.
2. Do not pipet by mouth, inhale, ingest or allow product to come into contact with open wounds. Wash thoroughly any part of the body which comes into contact with the product.
3. Avoid accidental autoinjection by exercising extreme care when handling in conjunction with any injection device.
4. This product is intended for research use by qualified personnel only. It is not intended for use in humans or as a diagnostic agent. Advanced Targeting Systems is not liable for any damages resulting from the misuse or handling of this product.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols



KNRK-NK1r cells were plated at 2500 cells/well and incubated overnight. Saporin (Cat. #PR-01), SSP-SAP (Cat. #IT-11) and Blank-SAP (Cat. #IT-21) were added in 10 ul volumes and the plates left to incubate for 72 hours. PMS/MTS developing reagent was added and the plates read at 490 nm.