

Targeting Trends

Reporting the latest news in Molecular Surgery

Cerebral cholinergic lesion reduces operant responses to unpleasant thermal stimuli

by Ronald G. Wiley, M.D., Ph.D., Departments of Neurology and Pharmacology, Vanderbilt University, Nashville, TN and C. J. Vierck, Department of Neuroscience, McKnight Brain Institute, University of Florida College of Medicine, Gainesville, FL, USA

Degeneration of the cholinergic basal forebrain (CBF: medial septum, diagonal band of Broca, nucleus basalis of Meynert/substantia innominata) is a prominent feature of Alzheimer's disease (AD). The CBF supplies cholinergic input to most of the cerebral cortex and hippocampus including somatosensory areas and anterior cingulate cortex that are involved in pain perception and experiencing discomfort, respectively. Clinical literature suggests that patients with AD either feel less pain or express discomfort less than comparable patients without dementia. As a result, AD patients receive less analgesics, but there is concern that AD only impairs communicating discomfort. Rats with extensive CBF lesions show impairment in a wide range of learning tasks and ability to sustain selective arousal/attention, but it is not known what role the CBF plays in central pain processing.

The present study sought to assess the impact of CBF lesions on behavioral responses to nociceptive stimuli in rats. Rats were trained on a thermal escape task where they chose whether to spend time in a dark chamber with the floor temperature at 10° C or 44.5° C (both mildly unpleasant), or move to a connected room temperature chamber with bright lighting. After establishing baseline performance on the operant task, selective CBF lesions were produced by intracerebroventricular injection of 192-saporin (192-IgG-SAP, Cat. #IT-01; Fig. 1). This immunotoxin selectively destroys neurons expressing p75NTR, the low affinity neurotrophin receptor that is uniquely expressed by CBF neurons in the forebrain. The rats were retested repeatedly over 19 weeks. On several occasions, the rats were subjected to sound stress prior to escape testing, and





Fig. 1: Representative coronal (frontal) sections from control rats (A, C) and 192-sap-treated rats (B, D) showing loss of CBF cholinergic neurons in the 192-sap rats. (A, B) Show the medial septal nucleus. (C, D) Show the nucleus basalis/substantia innominate regions. Sections were stained for demonstration of choline acetyltransferase using the immunoperoxidase technique (see text). The magnification bars in lower right corners indicate 250 μm (1).

(continued on page 6)

 Brian Russell - Editor	Letter from the President: "A sigh is (not) just a sigh..." Page 2	New Products & Tools: Vesicular GABA Transporter Transfected cell lines Page 3	Journal Time: Latest Pubs & Refs reviewed Pages 4-5	Talking About Targeting: Saporin Safety – Nothing to Fear Page 7
---	---	--	--	---

A Sigh is (Not) Just a Sigh . . .

Denise Higgins - President



*You must remember this
A kiss is just a kiss,
A sigh is just a sigh.
The fundamental things apply
As time goes by.*

Casablanca theme song, that prompted the famous line, "Play it again, Sam." If you haven't seen this classic 1944 Academy Award winner, I highly recommend it!

The fundamental things apply. "The Peptidergic Control Circuit for Sighing," recently published in the prestigious journal *Nature*, has made us rethink our fundamental belief that sighs are only "long, deep breaths expressing sadness, relief or exhaustion." Often prompting someone to say, "What's wrong?" As it turns out, sighs "also occur spontaneously every few minutes to reinflate alveoli, and sighing increases under hypoxia, stress, and certain psychiatric conditions." Thanks to the clever researchers led by Dr. Jack Feldman at UCLA, and their collaborators at Stanford University School of Medicine, we now know a lot more about this process (see Fig. 1 below and Reference Summary on Pg. 4).

The Bötzing Complex plays an important role in controlling breathing and was named by UCLA Professor Jack Feldman in 1978. Feldman named this area after a bottle of white wine named Botzinger present at his table (perhaps he was allowing it to breathe) during a scientific meeting in Hirschhorn, Germany, that year. Jack Feldman named the most rostral portion of the ventral respiratory group and continues to pave the way for important respiratory research. It's a song we all need to hear, so:

Play it again, . . . Jack!

Deservedly, Jack Feldman's findings went viral. Here are just a few of the links to interviews and commentary:

NEWS ARTICLES

The Washington Post: "Scientists uncover the brain mechanism that makes you sigh."

The LA Times: "Scientists locate the part of the brain where sighs are made."

NPR: "Sorry, Bogie, a Sigh Is Not Just A Sigh."

WSJ: "Scientists Pinpoint Brain Chemical Linked to the Sigh."

THE GUARDIAN: "A sigh's not just a sigh -- it's a fundamental life-sustaining reflex."

THE ONION: "Sighing a Life-Sustaining Reflex."

RADIO INTERVIEWS BBC: "How the brain's sighing reflex was named."

NZ Radio: "Sigh science."

TELEVISION COVERAGE CBS This Morning: "More than just a sigh."

Bötzing Complex References Using ATS Products

Gray PA, Janczewski WA, Mellen N, McCrimmon DR, Feldman JL. (2001) Normal breathing requires preBotzinger complex neurokinin-1 receptor-expressing neurons. *Nat Neurosci* 4(9):927-930 (SP-SAP, Cat. #IT-07).

Feldman JL, Mitchell GS, Nattie EE (2003) BREATHING: Rhythmicity, Plasticity, Chemosensitivity. *Annu Rev Neurosci* 26:239-266 (SERT-SAP, Cat. #IT-23; SP-SAP, Cat. #IT-07*).

Wenninger JM, Pan LG, Klum L, Leekley T, Bastastic J, Hodges MR, Feroah T, Davis S, Forster HV (2004) Small reduction of neurokinin-1 receptor-expressing neurons in the pre-Botzinger complex area induces abnormal breathing periods in awake goats. *J Appl Physiol* 97(5):1620-1628 (SP-SAP, Cat. #IT-07*).

Wenninger JM, Pan LG, Klum L, Leekley T, Bastastic J, Hodges MR, Feroah TR, Davis S, Forster HV (2004) Large lesions in the pre-Botzinger complex area eliminate eupneic respiratory rhythm in awake goats. *J Appl Physiol* 97(5):1629-1636 (SP-SAP, Cat. #IT-07*).

McKay LC, Janczewski WA, Feldman JL (2005) Sleep-disordered breathing after targeted ablation of preBotzinger complex neurons. *Nat Neurosci* 8(9):1142-1144 (SP-SAP, Cat. #IT-07*).

McKay LC, Feldman JL (2008) Unilateral Ablation of preBotzinger Complex Disrupts Breathing During Sleep but not Wakefulness. *Am J Respir Crit Care Med* 178(1):89-95 (SP-SAP, Cat. #IT-07*).

Montandon G, Qin W, Liu H, Ren J, Greer JJ, Horner RL. (2011) PreBotzinger complex neurokinin-1 receptor-expressing neurons mediate opioid-induced respiratory depression. *J Neurosci* 31(4):1292-1301 (anti-NK1r Cat. #AB-N04**).

Gray PA, Hayes JA, Ling GY, Llona I, Tupal S, Picardo MCD, Ross SE, Hirata T, Corbin JG, Eugenin J, Del Negro CA (2010) Developmental origin of preBotzinger Complex respiratory neurons. *J Neurosci* 30(44):14883-14895 (anti-NK1r Cat. #AB-N04**).

*See alternate product: SSP-SAP (Cat. #IT-11); **See alternate product: NK-1r affinity purified antibody (Cat. #AB-N33AP)




Fig. 1 On each side of the brain stem, a fluorescent-green marker illuminates the 200 neurons that control the sighing reflex. To determine if Neuromedin B receptor (NMBR)- and Gastric Releasing Peptide receptor (GRPR)-expressing neurons function specifically in sigh control, they were removed using Bombesin-SAP (Cat. #IT-40); Bombesin binds both receptors.

Photo Credit: Stanford/Krasnow lab

ATS Toolbox - New Products

vGAT Products

ATS is pleased to present a new product line specific for the vesicular GABA transporter (vGAT) protein. vGAT mediates both the accumulation of GABA into synaptic vesicles and its release from nerve terminals. vGAT is expressed in nerve endings of GABAergic neurons throughout the CNS. The GABAergic system is crucial to the development and functional maturation of the nervous system, as well as the maintenance of balance between excitation and inhibition required for normal neural circuit function. Dysfunction of GABAergic neurons underlies aspects of clinical symptoms found in several diseases such as epilepsy, Down Syndrome, Fragile X Syndrome, Schizophrenia, and Autism among others.

Anti-vGAT-SAP (Cat. #IT-71) is highly specific for cells that express vGAT. Instead of spending precious time and money producing a vGAT knockout animal, you can use Anti-vGAT-SAP to specifically eliminate cells that express vGAT. Anti-vGAT-SAP also allows you to study the behavioral effects before and after treatment and subsequent elimination of vGAT expressing cells. We have demonstrated the specificity of our Anti-vGAT-SAP by cytotoxicity assays (Fig. 1).




Fig. 1: Cells were plated at 1000 cells/90 μ l/well in a 96 well plate and incubated overnight. Anti-vGAT-SAP was added in 10 μ l volumes and the plates were incubated for 72 hours. The plates were developed with SRB and read at 564 nm in a plate reader. Data analysis was done by PRISM (GraphPad, San Diego).


In addition to the targeted toxin, we are offering other vGAT products for your ELISA, flow cytometry and immunoblotting needs. We currently offer rabbit polyclonal vGAT antibody serum (AB-N44), affinity-purified antibody (AB-N44AP), biotinylated antibody (BT-N44), and Alexa-488 conjugated antibody (FL-N44). Keep an eye out for our transfected vGAT cell line, coming soon!

Stable Transfected Cell Lines

Advanced Targeting Systems is proud to announce that we are now making available to the public our proprietary line of stably transfected cell lines, used for years in-house to validate ATS targeted toxins:

192 IgG-SAP was the first ATS product and lesioning use of the product has been published 400+ times over 20+ years, so it is only appropriate that we launch this new product line with CHO_{p75} cells, now available for purchase. Used for some time now as the preferred cell line for QC verification of each new lot of 192 IgG-SAP, CHO_{p75} cells are standard CHO (Chinese Hamster Ovary) cells that stably express the rat low affinity nerve growth factor, p75 (p75^{NTR}). Expression of p75^{NTR} in CHO_{p75} cells has been verified by flow cytometry with 192-IgG-Alexa488 (Cat. #FL-03) and cytotoxicity assay with 192-IgG-SAP (Cat. #IT-01).

ATS is on the verge of releasing several other transfected cell lines, including those expressing VGAT and the orexin-2 receptor. Transfected cell lines are valuable tools that can be used for investigating the function of the transfected molecules. They are also highly useful for screening research and therapeutic reagents that target the transfected gene product. Visit the ATS website for updates and pricing on this new product offering.



CHO and CHO_{p75} cells were plated at 1000 cells/90 μ l/well and incubated overnight. 192-IgG-SAP and saporin were added in 10 μ l volumes, and the plates incubated for 72 hours. The cells were fixed with 10% TCA, then stained with 0.4% sulfarhodamine B/1% acetic acid. The plates were read at 564 nm. Data analysis done by Prism software (GraphPad, San Diego).

Recent Publications & References

Reviewed by Matthew Kohls

The peptidergic control circuit for sighing.

Li P, Janczewski WA, Yackle K, Kam K, Pagliardini S, Krasnow MA, Feldman JL.

Nature. 530(7590):293-297, 2016.

Sighs are often associated with relief or sadness, but rodents sigh spontaneously dozens of times per hour. There are physiological benefits to sighing, including enhancement of gas exchange and preservation of lung integrity. The authors identify a peptidergic sigh control circuit in the retrotrapezoid nucleus/parafacial respiratory group of the mouse brain that projects to the pre-Bötzinger complex. Mice received bilateral 6.2-ng injections of Bombesin-SAP (Cat. #IT-40) into the pre-Bötzinger complex. Blank-SAP (Cat. #IT-21) was used as control. Elimination of the bombesin receptor-expressing neurons or inhibition of neuromedin B receptor-expressing neurons suppressed sighing. Interfering with the activity of both receptors abolished sigh activity while leaving normal breathing intact. The data suggest that overlapping peptidergic pathways are the core of a sigh control circuit. (See article on Pg 2).

Pain sensitivity following loss of cholinergic basal forebrain (CBF) neurons in the rat.

Vierck CJ, Yeziarski RP, Wiley RG. *Neuroscience* 319:23-34, 2016.

There is a large amount of research on the involvement of cholinergic mechanisms on spinal transmission of pain signals, indicating that cholinergic agonists can attenuate this kind of pain. In contrast, some studies have shown affective reactions to pain are suppressed by cholinergic antagonists. The authors investigated the disagreement between reflexive and affective reactions with a 4- μ g 192-IgG-SAP (Cat. #IT-01) injection into the left lateral ventricle of rats. Animals were tested in temperature escape and sound stress models. Lesioned rats displayed decreased escape from thermal stimulation, as well as loss of the normal hyperalgesic effect of sound



stress. Results indicate that the basal forebrain cholinergic system plays a role in central processing of pain. (See Pg 1.)

Ablation of mu opioid receptor-expressing GABA neurons in rostromedial tegmental nucleus increases ethanol consumption and regulates ethanol-related behaviors.

Fu R, Chen X, Zuo W, Li J, Kang S, Zhou LH, Siegel A, Bekker A, Ye JH. *Neuropharmacology* 2016.

In this work the authors investigated cellular mechanisms underlying the aversive effects of alcohol that limit its intake. Previous work has linked synaptic inhibition of dopamine neurons in the ventral tegmental area to this aversion. Rats conditioned to ingest ethanol received bilateral injections totaling 3 pmol of Dermorphin-SAP (Cat. #IT-12) into the rostromedial tegmental nucleus (RTMg). Blank-SAP (Cat. #IT-21) was used as a control. Lesioned animals displayed significantly increased preference for, and intake of ethanol, while showing no change in the desire for sucrose. The results indicate that mu opioid expressing GABAergic neurons in the RTMg are highly involved in the regulation of ethanol consumption.

Substituting mouse transcription factor Pou4f2 with a sea urchin orthologue restores retinal ganglion cell development

Mao C-A, Agca C, Mocko-Strand JA, Wang J, Ullrich-Lüter E, Pan P, Wang SW, Arnone MI, Frishman LJ, Klein WH. *Proc Royal Society London B* DOI: 10.1098/rspb.2015.2978

Although the regulatory genes for eye development are highly conserved, there is a vision is widely diversified between species. Little is known about how gene networks vary to produce the variety of structures and functions seen across organisms. The authors investigated photoreception in echinoderms, adult sea urchins. Urchins have no structures resembling vertebrate eyes, but recent work has demonstrated the presence of photoreceptor neurons. In this work the authors transferred the urchin version of a transcription factor involved in retinal ganglion cell development into mice lacking the mouse version of that gene. The urchin gene was able to restore function in the mouse, indicating the depth of conservation for eye development gene networks. Some of the immunohistochemical staining was done with anti-melanopsin (Cat. #AB-N39) at a 1:1000 dilution.

Current and Future Issues in the development of spinal agents for the management of pain.

Yaksh TL, Fisher C, Hockman T, Wiese A. *Curr Neuropharmacol* 2016 Mar 7 [Epub]

Although conscious pain experience is driven by signals mediated supraspinally, the more high intensity pain generated by strong stimuli, tissue injury, and nerve injury is encoded at the spinal dorsal horn level. The control of pain signals at the spinal dorsal horn level is a tempting target for targeted pain therapy. This review discusses the potential targets for pain therapeutics in the spinal dorsal horn, and some of the spinal agents used to modulate pain transmission through that location. The use of SSP-SAP (Cat. #IT-11) is mentioned as a neurokinin-1 targeted

(continued on page 5)

Recent Publications & References

(continued from page 4)

molecule that can block some pain transmission.

Reorganization of Motor Cortex by Vagus Nerve Stimulation Requires Cholinergic Innervation.

Hulsey DR, Hays SA, Khodaparast N, Ruiz A, Das P, Rennaker RLn, Kilgard MP.

Brain Stimul 9(2):174-181, 2016.

Recent work has suggested that vagus nerve stimulation (VNS) can enhance neuroplasticity, and coupled with other training can drive motor cortex reorganization. These findings highlight the potential of VNS to support recovery from neurological disease. Pretrained rats received bilateral injections totaling 3.75 µg of 192-IgG-SAP (Cat. #IT-01) into the nucleus basalis (NB). Mouse-IgG-SAP (Cat. #IT-18) was used as control. Control animals displayed a substantial increase in proximal limb representation, lesion of the NB prevented this increase. Motor performance was similar between lesion and control groups, indicating that the difference in representation was not due to altered limb function.

The effect of nucleus basalis magnocellularis deep brain stimulation on memory function in a rat model of dementia.

Lee JE, Jeong da U, Lee J, Chang WS, Chang JW.

BMC Neurol 16(1):6, 2016.

Deep brain stimulation (DBS) is the application of electrical impulses to specific parts of the brain for treating disorders such as Parkinson's disease, chronic pain, and obsessive-compulsive disorder. It has been theorized that stimulation of brain structures associated with memory can enhance cognitive function. The authors lesioned the basal forebrain of rats through bilateral injections totaling 5 µg of 192-IgG-SAP (Cat. #IT-01) into the lateral ventricle. Animals then received DBS to the nucleus basalis magnocellularis and were tested in a Morris water maze task. Results indicate that DBS has beneficial effects on consolidation and retrieval of visuospatial memory.

Neuroplasticity and Repair in Rodent Neurotoxic Models of Spinal Motoneuron Disease.

Gulino R.

Neural Plast 2016:2769735, 2016.

TDP-43 (Transactive response DNA-binding protein) is a highly conserved nuclear protein that binds both DNA and RNA. It has been found in cytoplasmic protein aggregates of patients with conditions such as amyotrophic lateral sclerosis and Alzheimer's disease. In this work the authors examine the role of TDP-43 in spinal cord plasticity. Mice received bilateral 3-µg injections of CTB-SAP (Cat. #IT-14) into the lateral and medial gastrocnemius muscles. The results indicate that motor performance is dependent on expression of synapsin-I, which in turn may be dependent on TDP-43.

Aminopeptidase N (APN/CD13) as a target molecule for scirrhous gastric cancer.

Nohara S, Kato K, Fujiwara D, Sakuragi N, Yanagihara K, Iwanuma Y, Kajiyama Y.

Clin Res Hepatol Gastroenterol 2016 Jan 13 [Epub ahead of print]

Scirrhous gastric cancer has the worst prognosis of gastric carcinoma, and treatment with standard cancer therapies has had minimal success. In this work the authors target CD13 as a marker for scirrhous gastric cancer. A gastric cancer cell line was challenged with a CD13 antibody coupled to Mab-ZAP (Cat. #IT-04) in an in vitro cytotoxicity assay. The anti-CD13 complex was more cytotoxic than an anti-EpCAM-immunotoxin. These data, combined with flow cytometry analysis and enzyme activity assays, demonstrate the expression of CD13 as a marker for scirrhous gastric cancer.

Perinatal 192 IgG-Saporin as Neuroteratogen.

Petrosini L, De Bartolo P, Cutuli D, Gelfo F. *Curr Top Behav Neurosci* 2015 Dec 23. [Epub ahead of print].

The authors discuss the effects of perinatal administration of 192-IgG-SAP (Cat. #IT-01) and areas of research that have been investigated through the use

of these lesions. The chapter covers a description of 192-IgG-SAP, lesioning methods, and outlines the short- and long-term biochemical, structural, behavioral, and cognitive effects of 192-IgG-SAP administration.

Neuroteratology and Animal Modeling of Brain Disorders.

Archer T, Kostrzewa RM.

Curr Top Behav Neurosci 2016 Feb 9. [Epub ahead of print].

This work covers development and use of the neurotoxins that are most commonly used as neuroteratologic agents - producing permanent, lifelong destruction of specific groups of neurons. Saporin conjugates are discussed, in terms of animal models of human neurodegenerative, neuropsychiatric, and neurological conditions.

Locus Coeruleus and Tubermammillary Nuclei Ablations Attenuate Hypocretin/Orexin Antagonist-Mediated REM Sleep.

Schwartz MD, Nguyen AT, Warriar DR, Palmerston JB, Thomas AM, Morairty SR, Neylan TC, Kilduff TS.

eNeuro. 2016 Mar 21;3(2). pii: ENEURO.0018-16.2016.

To examine the mechanism by which the Orexin 1r/Orexin 2r antagonist almorexant decreases wakefulness and increases NREM and REM sleep the authors utilized Anti-DBH-SAP (Cat. #IT-03) and Orexin-SAP (Cat. #BETA-031). Rats received 3-µg injections of Anti-DBH-SAP into the LC, or bilateral 57-80 ng injections of Orexin-SAP into the TMN. Both conjugates attenuated the increased REM sleep seen upon administration of almorexant without altering almorexant-induced changes in NREM sleep.

Did we miss you?
If so, please send a PDF of your publication to
admin@atsbio.com so we can feature your research in our next issue!

Cerebral cholinergic lesion reduces operant responses to unpleasant thermal stimuli

(continued from page 1)

the rats were also tested on the thermal plate (hot/cold plate) without an escape option to measure lick guard (reflex) responses with and without preceding stress.

Compared to controls, 192-IgG-SAP injected rats showed highly significant ($p < 0.001$) loss of neurons from all subdivisions of the CBF based on post mortem brain sections stained for choline acetyltransferase. The CBF-lesioned rats escaped less than controls after 192-IgG-SAP injection (i.e. less motivated to get away from the aversively hot or cold stimuli). Reflex lick/guard responses, which are mediated at the spinal level, were not affected. The usual hyperalgesic effect of stress on the operant thermal escape task was absent in the CBF-lesioned rats. These results indicate a role for the CBF in modulating central pain processing. The loss of stress effect on thermal escape responses is consistent with loss of the arousal/attention function(s) of the CBF. These results also demonstrate the usefulness of 192-IgG-SAP for studies of the role of central (cerebral) cholinergic mechanisms in pain processing and are consistent with the idea that AD patients experience less discomfort for a given painful condition.

Reference

1. Pain sensitivity following loss of cholinergic basal forebrain (CBF) neurons in the rat. Vierck CJ, Yeziarski RP, Wiley RG. *Neuroscience* 319:23-34, 2016.

Where to See Us:

Am. Assoc. of Cancer Research
New Orleans, April 17-20

Am. Assoc. of Immunology
Seattle, May 13-17


American Pain Society
Austin, May 11-14

Targeting Membrane Proteins
Boston, June 16-17

For years ATS fans have been Gangsta and Kermit fans, and now it is time to show your love by adding a caption to Kermit's picture! Each quarter we will post a picture of Kermit up to his usual hijinks and ask our readership to offer their most creative captions.



Winning Caption Goes Here!!!



Be clever and creative as Kermit is judging you from up on high!

Email your captions to:

admin@atsbio.com
Subject line:
Kermit Caption.

The winning caption will be posted in the next newsletter and the author will receive a

\$350 product credit
runner-up will receive
\$100 product credit

Talking about Targeting

Saporin Safety

Over the years, ATS has frequently been confronted with questions about Saporin's safety for use in the lab as well as when used clinically. Residual awareness of alternate Ribosome-Inactivating Proteins (RIPs) and 'toxins' such as Ricin have poisoned the belief that Saporin is safe. As a Type I RIP, Saporin has no binding chain and consequently no means of entering the physiological space necessary for the protein to act as a toxin. As such, in response to specific concerns about safety from casual users, reviewers of work with saporin, and potential 3rd-party manufacturers of saporin and SP-SAP, the following is a review of safety in handling and potential toxicity within the human body for systemic events not related to the predicted therapeutic application of SP-SAP.

The acute LD50 for saporin in mice (25 g) is 6.8 mg/kg;¹ that would translate in humans (75 kg) to 510 mg! A concentration of about 100 nM is the threshold to see even a vague hint of saporin toxicity. In human blood, that would correspond to 24 mg injected systemically into a person. The fermentation process to produce recombinant saporin has a titer of 2 mg/L meaning that the production broth itself contains no more than 67 nM concentration of saporin. Furthermore, the final protein concentrations from production batches of recombinant Saporin used in our drug are 4 mg/ml, meaning 6 mL of final material would need to accidentally end up in a human before the 'hint of toxicity' threshold would potentially be met.

The toxicology studies of SP-SAP contained within ATS's IND prior to the current human Phase I clinical trial evaluated effects related to the intended method of administration, intrathecal local injection. SP-SAP is not expected to ever be a self-administered therapy, so the effects of gross off-target events, such as accidental auto-injection, swallowing, spillage, or immersion were not considered.

The table below² highlights antibody-saporin conjugates approved by the FDA for Phase I/II clinical trials in humans. The therapeutics listed below were administered intravenously and imply what the FDA accepted as non-toxic levels of saporin-based conjugates in these studies.

Table 2. Clinical trials in patients with SAP containing ITs.

Antibody	Antigen	Disease	Total Dose	PR	SD/MR	No. patients	Ref.
F(ab') ₂ BsAb	CD22	NHL	5 mg	-	1	1	[32]
4KB128 + HD6	CD22	B-cell lymphoma	5-20 mg	-	4	4	[33]
F(ab') ₂ BsAb	CD22	NHL	5-20 mg	-	5	5	[34]
Ber-H2	CD30	HD	0.8 mg/kg	3 (75%)	1 (25%)	4	[35]
Ber-H2	CD30	HD	0.2-0.8 mg/kg	5 (40%)	3 (25%)	12	[36]

HD: Hodgkin's disease, NHL: non-Hodgkin's lymphoma, PR: partial remission, SD/MR: stable disease/minor response.

Looking more closely at the study by French *et al.*,³ several milligrams of antibody conjugate were repeatedly injected into human patients under a FDA regulated clinical trial and peak serum levels tested, demonstrating rapid clearing of saporin from the system.

As a company that specializes in Saporin, our two-plus decades of experience working with the protein in research, preclinical, and clinical environments has taught us that with minimal standard laboratory precautions users are not at any real risk of toxic effects. Even our CSO, after 30+ years of working with Saporin exhibits undetectable levels of Saporin antibodies in his blood!

References

1. Thorpe, P.E.; Brown, A.N.; Bremner, J.A., Jr.; Foxwell, B.M.; Stirpe, F. An immunotoxin composed of monoclonal anti-Thy 1.1 antibody and a ribosome-inactivating protein from *Saponaria officinalis*: Potent antitumor effects *in vitro* and *in vivo*. *J. Natl. Cancer Inst.* 1985, 75, 151-159.
2. Polito L, Bortolotti M, Pedrazzi M, Bolognesi A. Immunotoxins and other conjugates containing saporin-s6 for cancer therapy. *Toxins* (Basel). 2011 Jun;3(6):697-720.
3. French, R.R.; Bell, A.J.; Hamblin, T.J.; Tutt, A.L.; Glennie, M.J. Response of B-cell lymphoma to a combination of bispecific antibodies and SAP. *Leuk. Res.* 1996, 20, 607-617.

Tollfree Phone: (877) 889-2288
Phone: (858) 642-1988


www.ATSbio.com

Monday-Friday 9am-5pm (PST).
Order online, via telephone or
fax, or with your local distributor.



10451 Roselle Street #300
San Diego, CA 92121


PRESORTED STD.
U.S. POSTAGE
PAID
SAN DIEGO, CA
PERMIT # 2686



A targeting agent is conjugated to a payload. here the conjugate is an antibody linked to a small molecule.



The targeted conjugate binds only the cells that are specific to the antibody causing the receptor to internalize.



The internalized conjugate is digested by enzymes releasing the payload into the cytosol. The payload is then free to alter the cell.

Daily Promotions on Targeting Tools

• Check our website daily •

Specials from each of our Product Managers



Kit Upgrade

Purchase a 250-mcg size Targeted Conjugate and receive a FREE upgrade to a KIT that includes all the specific controls necessary to validate your assay.

Enter **TTKIT-UP** at checkout. Offer expires June 30, 2016.

Brian Russell

Free Antibody with Transfected Cell Line

CHOp75 cells stably express the rat low affinity nerve growth factor, p75 (p75NTR). Order this cell line and receive 1 mg of positive control antibody **FREE** until June 30, 2016.



Matt Kohls



Save on ZAP products

Let ATS help you screen your antibodies with our secondary conjugate product line. Choose your size: 100-mcg for 75-mcg price; 400-mcg for 250-mcg price; 1000-mcg for 500-mcg price. Use Discount Code **ZAPQ216**.

Offer expires June 30, 2016.

Leonardo Ancheta

Biotin or Fluorescent labeling

Order one Custom labeling and get a second one for 50% off! OR Buy 2 services and get the 3rd service FREE! With the purchase of a biotinylation service we even include a free Biotin-Z Kit to verify internalization of your biotinylated material. Offer expires June 30, 2016. Use promo code **CustomBioFluo**.



Patrick Shramm



Dopamine makes you feel good . . .

but so does saving lots of money! Order 50 mcg of Anti-DBH-SAP, Anti-DAT-SAP or the corresponding kits and double your order to 100 micrograms for free! But it doesn't end there. Purchase 2 vials of a dopamine-related antibody or antigen and get a third for free! Enter coupon code **DOPAX2** through Offer expires June 30, 2016.

Chelsea Friedman