

**Anti-Human CXCR4-SAP**
TARGETED SAP CONJUGATE*[affinity-purified anti-human CXCR4]-streptavidin-saporin*

Catalog Number: BETA-018
Quantity: 25 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative. Sterile-filtered.

Background: Targeted SAP conjugates are powerful and specific lesioning agents used in the technique known as Molecular Surgery. The ribosome-inactivating protein, saporin (from the seeds of the plant, *Saponaria officinalis*) is bound to a targeting agent (anything that is recognized on the cell surface and internalized). The targeted conjugate is administered to cells (*in vitro* or *in vivo*). The targeting agent seeks out and binds to its target on the cell surface. The conjugate is internalized, saporin breaks away from the targeting agent, and inactivates the ribosomes which causes protein inhibition and, ultimately, cell death. Cells that do not have the cell surface marker are not affected.

The chemokine receptor type 4 (CXCR4, fusin, CD184), is a G-protein-coupled receptor consisting of a 7-transmembrane region and is the primary receptor for stromal cell-derived factor-1 (SDF-1) also called CXCL12. CXCR4 has become a topic of interest with its involvement in many different cancer types such as breast and lung in particular with implications in ovarian, prostate, colon, ovary and bladder as well. CXCR4 has also been characterized in its involvement as the main receptor supporting viral entry for HIV and evidence also showing that CXCR4 mutations causes WHIM Syndrome, a group of immunodeficiency disease which include warts, Hypogammaglobulinemia, infections and Myelokathexis.

Specificity & Preparation: This targeted toxin recognizes cells that express human chemokine receptor type 4 (CXCR4, fusin, CD184). Anti-human CXCR4-SAP is a bonded conjugate of a mouse monoclonal antibody to human CXCR4 and the secondary conjugate Streptavidin-ZAP containing the ribosome-inactivating protein, saporin.

Usage: Anti-human CXCR4-SAP eliminates cells expressing human CXCR4 (fusin, CD184). All other cells are left untouched. It is useful in retrograde transport (see Wiley *et al*, 1989). **There may be lot-to-lot variation in material; working dilutions must be determined by end user. If this is a new lot, you must assess the proper working dilution before beginning a full experimental protocol.**

Storage: Gently spin down material 5-10 seconds in a microfuge before use. Store the material in undiluted aliquots at -20°C . Material should be aliquoted to a convenient volume and quantity to avoid repeated freezing and thawing that can damage the protein content. Under these conditions, the material has a very stable shelf-life. Thawing should be done at room temperature or on ice. The thawed solution should remain on ice until use.

Do not use a reducing agent (such as dithiothreitol, beta-mercaptoethanol or ascorbic acid) with this material. It will inactivate the toxin.

This material is an extremely potent cytotoxin. Handling should be done by experienced personnel. Gloves and safety glasses are required when handling this product. Care in disposal is mandatory; autoclaving or exposure to 0.2 M sodium hydroxide will inactivate the material. All labware that comes into contact with this material should be likewise treated.

**Anti-Human CXCR4-SAP**
TARGETED SAP CONJUGATE**Selected References:**

1. Wiley RG, Stirpe F, Thorpe P, Oeltmann TN (1989) Neuronotoxic effects of monoclonal anti-Thy 1 antibody (OX7) coupled to the ribosome inactivating protein, saporin, as studied by suicide transport experiments in the rat. *Brain Res* 505:44-54.

Control(s): BIgG-SAP Mouse**Safety:**

Good laboratory technique must be employed for safe handling of this product. This requires observation of the following practices:

1. Wear appropriate laboratory attire, including lab coat, gloves and safety glasses.
2. Do not pipet by mouth, inhale, ingest or allow product to come into contact with open wounds. Wash thoroughly any part of the body which comes into contact with the product.
3. Avoid accidental autoinjection by exercising extreme care when handling in conjunction with any injection device.
4. This product is intended for research use by qualified personnel only. It is not intended for use in humans or as a diagnostic agent. Advanced Targeting Systems is not liable for any damages resulting from the misuse or handling of this product.