



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to Acetylated Lysine
Rabbit Polyclonal

Catalog Number: AB-202
Quantity: 100 micrograms
Format: TBS, 0.1% NaN₃
Host: Rabbit
Immunogen: acetylated KLH conjugates

Specificity and Preparation:

Recognizes proteins acetylated on lysine residues. Tested: acetylated histone, acetylated BSA, and acetylated MBP, no reaction to the non-acetylated proteins. Purified with immobilized acetylated lysine on agarose. Concentration 250 μ g/ml.

Usage and Storage:

Applications include ELISA, Western Blot, and immunoprecipitation. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.