



11175-A FLINTKOTE AVENUE, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Antibody to AKT2 (PKB-beta)  
Rabbit Polyclonal**

**Catalog Number:** AB-205  
**Quantity:** 100 micrograms  
**Format:** 250  $\mu\text{g}/\text{mL}$  in TBS pH 7, 0.01% NaN<sub>3</sub>  
**Host:** Rabbit  
**Immunogen:** peptide(ASIRE) based on human PKB beta-KLH conjugates

**Specificity and Preparation:**

Recognizes PKB beta, does not cross-react with PKB alpha and gamma. Immunoaffinity purified with PKB beta peptide (ASIRE) on Agarose.

**Usage and Storage:**

Applications include ELISA (0.25  $\mu\text{g}/\text{mL}$ ) and Western Blot (1:1000). Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.