



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to AKT2 (PKB-beta)
Rabbit Polyclonal

Catalog Number: AB-205
Quantity: 100 micrograms
Format: TBS, pH 7, 0.01% NaN₃
Host: Rabbit
Immunogen: peptide(ASIRE) based on human PKB beta-KLH conjugates

Specificity and Preparation:

Recognizes PKB beta, does not cross-react with PKB alpha and gamma. Immunoaffinity purified with PKB beta peptide (ASIRE) conjugated to agarose. The concentration is 250 µg/ml.

Usage and Storage:

Reported to be effective for ELISA (0.25 µg/mL) and immunoblotting (western blot, 1:1000).

Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols