



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Anti-Phospho-S6 Peptide
Rabbit Polyclonal

Catalog Number: AB-207
Quantity: 100 micrograms
Format: PBS, pH 7, 0.01% NaN₃
Host: Rabbit
Immunogen: phosphorylated S6 peptide KRRRLApSLR-KLH conjugates

Specificity and Preparation:

Recognizes only the phosphorylated S6 peptide product of S6K, does not cross-react with the non-phosphorylated S6 peptide substrate. Affinity-purified with phospho-peptide KRRRLApSLR on Agarose and affinity-adsorbed with non-phosphorylated S6 peptide. Concentration 250 µg/ml.

Usage and Storage:

Applications include Western Blot, ELISA, and non-radioactive S6K kinase assay. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.