



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to HA-tag
Mouse Monoclonal

Catalog Number: AB-212
Quantity: 100 micrograms
Format: PBS, pH 7.4 with 0.05% sodium azide. Frozen liquid.
Host: Mouse
Isotype: IgG3
Clone: HA.C5
Immunogen: A synthetic peptide from influenza hemagglutinin epitope (YPYDVPDYA) coupled to KLH

Specificity and Preparation:

Recognizes HA-tagged proteins over-expressed in cells, including both amino- or carboxy-termini of targeted proteins in transfected mammalian cells. Concentration 1 mg/ml.

Usage and Storage:

Applications include Western Blot (1:1000), immunoprecipitation and immunostaining (1:200). Store the antibody at 4°C for three months or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.