



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Antibody to Histone H4**  
Rabbit Polyclonal

**Catalog Number:** AB-216  
**Quantity:** 100 micrograms  
**Format:** TBS, pH 7.0, 0.01% NaN<sub>3</sub>  
**Host:** Rabbit  
**Immunogen:** human histone H4 c-terminal peptide KLH conjugates

**Specificity and Preparation:**

Recognizes histone H4 in human, rat, bovine, mouse, and rabbit. Does not cross-react with histone H1, H2 and H3. Affinity-purified with histone H4 peptide conjugated to agarose. The concentration is 250 µg/ml.

**Usage and Storage:**

Reported to be effective for ELISA, immunoblotting (western blot), and immunoprecipitation. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

**To view protocol(s) for this and other products please visit: [www.ATSBio.com/protocols](http://www.ATSBio.com/protocols)**