



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to Human Serum Albumin (HSA)
Rabbit Polyclonal

Catalog Number: AB-218
Quantity: 1 milligram
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: human serum albumin whole molecule (HSA)

Specificity and Preparation:

Recognizes 70 kD human albumin, cross-reacts with bovine serum albumin, but not with human IgG. Affinity-purified with HSA on Agarose followed by affinity adsorption with human IgG. Concentration 2.5 mg/ml.

Usage and Storage:

Applications include ELISA, Western Blot, immunoprecipitation, and urinary protein capture assay. Store the antibody at -20°C. Stable for up to 36 months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.