



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

HRP-labeled Antibody to JNK
Rabbit Polyclonal

Catalog Number: AB-219
Quantity: 100 micrograms
Format: Stabilizer with 50% glycerol
Host: Rabbit
Immunogen: synthetic peptide based on human JNK1 (175-189 aa)-KLH conjugates

Specificity and Preparation:

Recognizes JNK. The recognition sequence is identical in JNK3, JNK2, and JNK1, in human, mouse, and rat. Affinity-purified with peptide-affinity chromatography. Conjugated to horse radish peroxidase (HRP) via reductive amination. Concentration 250 $\mu\text{g/ml}$.

Usage and Storage:

Applications include ELISA (1:1000) and Western Blot (1:250-1:1000). Store the antibody at -20°C . Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.