



11175-A FLINTKOTE AVENUE, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Antibody to JNK
Rabbit Polyclonal**

Catalog Number: AB-220
Quantity: 100 micrograms
Format: 250 μ g/mL in stabilizer with 50% glycerol
Host: Rabbit
Immunogen: synthetic peptide based on human JNK1 (175-189 aa)-KLH conjugates

Specificity and Preparation:

Recognizes JNK. The recognition sequence is identical in JNK3, JNK2, and JNK1, in human, mouse, and rat. Affinity-purified with peptide-affinity chromatography. Conjugated to horse radish peroxidase (HRP) via reductive amination.

Usage and Storage:

Applications include ELISA (1:1000) and Western Blot (1:250-1:500). Store the antibody at -20°C. Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.