



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

HRP-labeled Anti-p33 ING1
Rabbit Polyclonal

Catalog Number: AB-229
Quantity: 100 micrograms
Format: TBS, 0.1% NaN₃
Host: Rabbit
Immunogen: mouse p33 ING1 N-terminal peptide-KLH conjugates

Specificity and Preparation:

Recognizes p33 ING1 of human and mouse. Does not cross-react with other ING family. Affinity-purified with the p33 ING c-terminal peptide on Agarose. Conjugated to horse radish peroxidase (HRP) via reductive amination. Concentration 250 $\mu\text{g/ml}$.

Usage and Storage:

Applications include ELISA (0.25 $\mu\text{g/mL}$) and Western Blot (1 $\mu\text{g/mL}$). Store the antibody at 4°C. Stable for up to six months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.