



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Antibody to p53 (ser-37)**  
Rabbit Polyclonal

**Catalog Number:** AB-238  
**Quantity:** 100 micrograms  
**Format:** PBS, pH 7 with 50% glycerol  
**Host:** Rabbit  
**Immunogen:** p53 peptide (ser-37)-KLH conjugates

**Specificity and Preparation:**

Recognizes only phosphor p53 in human, not mouse. Affinity-purified with peptide-affinity chromatography. Concentration 250  $\mu\text{g/ml}$ .

**Usage and Storage:**

Applications include ELISA (1:2500) and Western Blot (1:2500). Store the antibody at  $-20^{\circ}\text{C}$ . Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.