



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Biotin-labeled Anti-p70 S6 Kinase**  
Rabbit Polyclonal

**Catalog Number:** AB-239  
**Quantity:** 100 micrograms  
**Format:** TBS, 0.1% NaN<sub>3</sub>  
**Host:** Rabbit  
**Immunogen:** human p70 S6K c-terminal peptide (ISKRPEHLRMNL)-KLH conjugates

**Specificity and Preparation:**

Recognizes p70 S6K of human, mouse, and rat. Does not cross react with p90 RSK. Affinity-purified with peptide (ISKRPEHLRMNL) on Agarose. Conjugated to d-Biotin via acylation of primary amine with active ester of biotin (NHS-activated d-biotin). Concentration 250 µg/ml.

**Usage and Storage:**

Applications include ELISA and Western Blot. Store the antibody at 4°C. Stable for up to six months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.