



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Biotin-labeled Anti-Phospho-Crebtide**  
Rabbit Polyclonal

**Catalog Number:** AB-245  
**Quantity:** 100 micrograms  
**Format:** PBS with 50% glycerol  
**Host:** Rabbit  
**Immunogen:** phospho-Crebtide (KRREILSRRPpSYR)-KLH

**Specificity and Preparation:**

Recognizes only the phosphorylated Crebtide product of PKA and PKC in a non-radioactive kinase assay, does not cross-react with the non-phosphorylated Crebtide. Affinity-purified with KRREILSRRPpSYR on Agarose followed by affinity adsorption with non-phosphorylated Crebtide. Biotinylated with NHS-activated biotin and G25 purification. Concentration 250  $\mu\text{g/ml}$ .

**Usage and Storage:**

Applications include ELISA, Western Blot, and non-radioactive PKA and PKC kinase assay. Store the antibody at  $-20^{\circ}\text{C}$ . Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.