



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Anti-Phospho-Crosstide**  
Rabbit Polyclonal

**Catalog Number:** AB-249  
**Quantity:** 100 micrograms  
**Format:** PBS, pH 7, 0.01% NaN<sub>3</sub>  
**Host:** Rabbit  
**Immunogen:** phospho-Crosstide (PRAApTF)-KLH

**Specificity and Preparation:**

Recognizes only the phosphorylated Crosstide, the product of SGK, PKB in a non-radioactive kinase assay, does not cross-react with the non-phosphorylated Crosstide. Affinity-purified with PRAApTF on Agarose followed by affinity adsorption with non-phosphorylated Crosstide. Concentration 250  $\mu$ g/ml.

**Usage and Storage:**

Applications include ELISA, Western Blot, and non-radioactive AKT and SGK kinase assay. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.