



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Biotin-labeled Anti-Phosphotyrosine
Rabbit Polyclonal

Catalog Number: AB-256
Quantity: 100 micrograms
Format: PBS with 0.01% NaN₃
Host: Rabbit
Immunogen: iodoacetylphosphotyramine-KLH conjugates

Specificity and Preparation:

Pan-specific for phosphotyrosine, does not cross-react with phospho-serine or -threonine. This antibody has been used on both mouse and human samples. Affinity-purified with phosphotyrosine immobilized conjugated to agarose. Conjugated to d-biotin via acylation of the primary amine with active ester of biotin (NHS-activated d-biotin). The concentration is 250 µg/ml.

Usage and Storage:

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:1000). Store the antibody at 4°C. Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols