

Anti-Phospho-Tyrosine
Rabbit Polyclonal

Catalog Number: AB-258
Quantity: 100 micrograms
Format: PBS with 0.01% NaN₃
Host: Rabbit
Immunogen: iodoacetylphosphotyramine-KLH conjugates


Specificity and Preparation:

Recognizes proteins and peptides phosphorylated on tyrosine residues. Does not cross-react with phosphoserine and threonine. Affinity-purified with phosphotyrosine on agarose. The concentration is 250 μ g/ml.

Usage and Storage:

Reported to be effective for ELISA (1:1000), immunoblotting (western blot, 1:250), immunoprecipitation (10 μ g/500 μ g protein sample), immunohistochemistry (1:50), and non-radioactive tyrosine kinase assay. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols



Immunohistochemistry analysis of human breast cancer tissue with Anti-Phospho-Tyrosine.