



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to SCD-1
Mouse Monoclonal

Catalog Number: AB-259
Quantity: 100 micrograms
Format: PBS, pH 7.4 with 0.05% sodium azide. Frozen liquid.
Host: Mouse
Isotype: IgG2b
Clone: CD.E10
Immunogen: full-length recombinant protein

Specificity and Preparation:

Reacts with human and mouse proteins. Concentration 1 mg/ml.

Usage and Storage:

Applications include Western Blot, ELISA, immunoprecipitation, and immunohistochemistry. Store the antibody at 4°C for three months or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.