



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Antibody to Cytochrome P450 (4F11)
MOUSE MONOCLONAL**

Catalog Number: AB-26
Quantity: 1 milliliter
Format: Hybridoma supernatant
Host: Mouse
Isotype: IgG1, k
Clone: F21 P6 F5
Immunogen: ovalbumin-conjugated synthetic peptide RVEPLGANSQ (c-terminal sequence)

Specificity and Preparation:

This antibody targets CYP P450 4F11 in homo sapiens. The target antigen is cytochrome P450 4F11 (EC 1.14.14.1) (CYP11B1). The immunogen was ovalbumin-conjugated synthetic peptide RVEPLGANSQ (c-terminal sequence) in BALB-c mouse.

Usage and Storage:

Applications include ELISA (undiluted, titer 1/1000), immunoblotting (western 1/10 dilution), and immunohistochemistry (1/5 dilution, antigen retrieval: microwave 20 min @ 800W in 10 mM citrate buffer, pH 6.0). Working dilutions must be determined by end user. Store the antibody at 4°C for one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

References:

1. Kumarakulasingham M, Rooney PH, Dundas SR, Telfer C, Melvin WT, Curran S, Murray GI (2005) Cytochrome p450 profile of colorectal cancer: identification of markers of prognosis. Clin Cancer Res 11 (10):3758-3765.