



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

FITC-labeled Rabbit Anti-Goat IgG
RABBIT POLYCLONAL

Catalog Number: AB-324
Quantity: 1 milligram
Format: 1 mg/mL in PBS, 50% Glycerol
Host: Rabbit
Immunogen: affinity-purified goat IgG, whole molecules

Specificity and Preparation:

Affinity-purified goat IgG antibodies are raised from rabbit using goat IgG whole molecules as an immunogen. The specific antibody is affinity purified with goat IgG immobilized on agarose and adsorbed with IgG from human and rabbit. The purified antibody is conjugated with FITC and the free FITC removed with G25 sephadex. Specificity is to goat IgG.

Usage and Storage:

The antibody is useful for detection, identification, quantitation, isolation and localization of goat IgG. Applications include ELISA (1:30,000), immunoblotting (western blot, 1:30,000), and immunohistochemistry (1:500).

Store at -20°C. Product is stable for up to 1 year.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols