



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**HRP-conjugated Antibody to HA
MOUSE MONOCLONAL**

Catalog Number: AB-466
Quantity: 500 microliters
Format: 1×PBS and 50% glycerol
Host: Mouse
Isotype: IgG2b
Clone: PHAHRPSHG
Immunogen: synthetic peptide (YPYDVDPDYA) coupled to KLH

Specificity and Preparation:

Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (YPYDVDPDYA) coupled to KLH.

Usage and Storage:

Applications include western blotting (1:500).

Antibody is stable for two years at -20°C, one month at 4°C.