



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Antibody to HA
RABBIT POLYCLONAL

Catalog Number: AB-469
Quantity: 100 microliters, 200 microliters, 500 microliters
Format: Serum with 0.02% sodium azide
Host: Rabbit
Clone: PPHASHG
Immunogen: synthetic peptide (YPYDVDPDYA) coupled to KLH

Specificity and Preparation:

Polyclonal antibodies are produced by immunizing rabbits with synthetic peptide (YPYDVDPDYA) coupled to KLH.

Usage and Storage:

Applications include western blotting (1:1,000) and immunoprecipitation.
Antibody is stable for two years at -20°C, one month at 4°C.