



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Anti-Conjugated Formiate
RAT POLYCLONAL**

Catalog Number: AB-T100
Quantity: 50 microliters
Format: Lyophilized and reconstituted with deionized water / 50% glycerol
Host: Rat
Isotype: IgG
Immunogen: Synthetic Formiate conjugated to bovine serum albumin

Specificity and Preparation:

Antiserum previously preabsorbed on protein carriers and purified by ammonium sulfate precipitation. This antibody targets conjugated Formiate. **This antibody does not recognize free formiate.**

Using a conjugate Formiate-BSA, antibody specificity was performed with an ELISA test by competition experiments with the following compounds:

COMPOUND	CROSS REACTIVITY §
Formiate-c-BSA	1
Lactate-c-BSA	1/>50,000
Butyrate-c-BSA	1/>50,000
Succinate-c-BSA	1/>50,000
Acetate-c-BSA	1/>50,000

BSA = Bovine Serum Albumin

§ Formiate-c-BSA concentration/conjugated close structure concentration at half displacement.

Usage and Storage:

Applications include ELISA (1/5,000-1/10,000). Store the antibody at 4°C for one month or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols