



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Conjugated NO-L-Methionine
ANTIGEN
NO-L-Methionine-G-BSA

Catalog Number: AG-079
Quantity: 250 micrograms, 1 milligram
Format: Lyophilized (1 mg); Lyophilized and reconstituted in deionized water (250 μ g)

Specificity and Preparation:

NO-L-Methionine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA), NO = nitrosylated.

Usage and Storage:

This antigen was used to produce a polyclonal antibody (Cat. #AB-T119). Applications include immunohistochemistry and immunocytochemistry.

Store at -20°C for one year. Reconstitute with deionized H₂O + 0.1% merthiolate (optional preservative).

Gently spin down material before use; 5-10 seconds in a microfuge should be adequate. This solution is stable at +4°C in darkness for 1 month.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols