



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121  
TELEPHONE (858) 642-1988 • FAX (858) 642-1989  
WWW.ATSBIO.COM • ATS@ATSBIO.COM

**Conjugated NO-Tyramine  
ANTIGEN**  
*NO-Tyramine-G-BSA*

**Catalog Number:** AG-080  
**Quantity:** 250 micrograms, 1 milligram  
**Format:** Lyophilized (1 mg); Lyophilized and reconstituted in deionized water (250  $\mu$ g)

**Specificity and Preparation:**

NO-Tyramine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA), NO = nitrosylated.

**Usage and Storage:**

This antigen was used to produce a polyclonal antibody (Cat. #AB-T123). Applications include immunohistochemistry and immunocytochemistry.

Store at -20°C for one year. Reconstitute with deionized H<sub>2</sub>O + 0.1% merthiolate (optional preservative). Gently spin down material before use; 5-10 seconds in a microfuge should be adequate. This solution is stable at +4°C in darkness for 1 month.

**To view protocol(s) for this and other products please visit: [www.ATSBio.com/protocols](http://www.ATSBio.com/protocols)**