



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

Conjugated Anthranilic Acid
ANTIGEN
Anthranilic Acid-BSA

Catalog Number: AG-121
Quantity: 250 micrograms, 1 milligram
Format: Lyophilized (1 mg); Lyophilized and reconstituted in deionized water (250 μ g)

Specificity and Preparation:
Anthranilic Acid conjugated with bovine serum albumin (BSA).

Usage and Storage:
This antigen was used to produce a polyclonal antibody (Cat. #AB-T155). Applications include immunohistochemistry and immunocytochemistry.
Store at -20°C for one year. Reconstitute with deionized H_2O + 0.1% merthiolate (optional preservative). Gently spin down material before use; 5-10 seconds in a microfuge should be adequate. This solution is stable at $+4^{\circ}\text{C}$ for 2 months.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols