



10451 ROSELLE STREET, #300, SAN DIEGO, CA 92121
TELEPHONE (858) 642-1988 • FAX (858) 642-1989
WWW.ATSBIO.COM • ATS@ATSBIO.COM

HRP-labeled Avidin

Catalog Number: PR-10
Quantity: 1 milligram
Format: 2 mg/mL in PBS buffer with 0.1% BSA

Specificity and Preparation:

Reacts with biotinylated substance. Conjugated to horseradish peroxidase (HRP) via reductive amination.

Usage and Storage:

Reported to be effective for ELISA, immunoblotting (western blot), and histochemistry.

Store the protein at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/protocols