

**Antibody to Acetylated Lysine**
RABBIT POLYCLONAL

Catalog Number: AB-202
Quantity: 100 micrograms
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: acetylated KLH conjugates

Specificity & Preparation: Recognizes proteins acetylated on lysine residues. Tested: acetylated histone, acetylated BSA, and acetylated MBP, no reaction to the non-acetylated proteins. Purified with immobilized acetylated lysine on agarose. The concentration is 250 µg/ml.

Usage: Reported to be effective for ELISA, immunoblotting (western blot), and immunoprecipitation.

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols