

HRP-labeled Anti-Histone H3
Rabbit Polyclonal

Catalog Number: AB-214
Quantity: 100 micrograms
Format: Stabilizing buffer, 50% glycerol
Host: Rabbit
Immunogen: human Histone H3 peptide (1-12)-KLH conjugates


Specificity and Preparation:

This antibody recognizes histone H3 in human and mouse. Affinity-purified with histone H3 peptide. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 250 mg/ml.

Usage and Storage:

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:250-1000). Store the antibody at -20°C. Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols



Western-blot analysis of Histone H3 expression in mouse spleen (Left) and mouse brain (right)