



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

HRP-labeled Antibody to JNK
Rabbit Polyclonal

Catalog Number: AB-219
Quantity: 100 micrograms
Format: Stabilizer with 50% glycerol
Host: Rabbit
Immunogen: Synthetic peptide based on human JNK1 (175-189 aa)-KLH conjugates

Specificity and Preparation:

Recognizes JNK. The recognition sequence is identical in JNK3, JNK2, and JNK1 in human, mouse, and rat. Affinity-purified with peptide-affinity chromatography. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 250 $\mu\text{g/ml}$.

Usage and Storage:

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:250-1:1000). Store the antibody at -20°C . Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols