



6352 CORTE DEL ABETO, STE B  
CARLSBAD, CA 92011 USA  
01.858.642.1988 • WWW.ATSBIO.COM

**Antibody to JNK**  
Rabbit Polyclonal

**Catalog Number:** AB-220  
**Quantity:** 100 micrograms  
**Format:** Stabilizer with 50% glycerol  
**Host:** Rabbit  
**Immunogen:** Synthetic peptide based on human JNK1 (175-189 aa)-KLH conjugates

**Specificity and Preparation:**

Recognizes JNK. The recognition sequence is identical in JNK3, JNK2, and JNK1, in human, mouse, and rat. Affinity-purified with peptide-affinity chromatography. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 250  $\mu\text{g/ml}$ .

**Usage and Storage:**

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:250-1:500). Store the antibody at  $-20^{\circ}\text{C}$ . Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

**To view protocol(s) for this and other products please visit: [www.ATSBio.com/support/protocols](http://www.ATSBio.com/support/protocols)**