



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Antibody to MBP (Myelin Basic Protein)
Rabbit Polyclonal

Catalog Number: AB-222
Quantity: 100 micrograms
Format: 0.1 M Tris-acetate, pH 7.0, 0.1% NaN₃
Host: Rabbit
Immunogen: bovine myelin basic protein (MBP)

Specificity and Preparation:

Recognizes 20 kD bovine myelin basic protein (MBP). May cross-react with other isoforms of MBP as the affinity binding-sequence is identical. Affinity-purified with pan MBP peptide (PRTPPPSQGKGRGL) conjugated to agarose. The concentration is 250 µg/ml.

Usage and Storage:

Reported to be effective for ELISA, immunoblotting (western blot), and immunoprecipitation. Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols