



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Biotin-labeled Anti-Methylated Lysine
Rabbit Polyclonal

Catalog Number: AB-223
Quantity: 100 micrograms
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: methylated KLH conjugates

Specificity and Preparation:

Recognizes proteins methylated on lysine residues (mono and di-methyllysine). Does not cross-react with acetylated proteins. Cross-reactivity with tri-methyllysine is not tested. Affinity-purified using Ne-monomethyllysine and dimethyllysine conjugated to agarose as affinity matrix. Conjugated to d-biotin via acylation. The concentration is 250 $\mu\text{g/ml}$.

Usage and Storage:

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:1000). Store the antibody at -20°C . Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols