



Antibody to Ochratoxin A RABBIT POLYCLONAL


Catalog Number: AB-227
Quantity: 100 micrograms
Format: PBS with 0.01% NaN₃
Host: Rabbit
Immunogen: ochratoxin A-KLH conjugates

Specificity & Preparation: Reacts with ochratoxin A and C, slightly cross-reacts with ochratoxin B, and does not cross-react with ochratoxin. Affinity-purified with ochratoxin A on agarose. The concentration is 250 µg/ml.

Usage: Reported to be effective for ELISA.

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols



Direct Competitive ELISA for Ochratoxin A
Immobilized anti-OA: 2 µg/mL
OA-HRP: 25 ng/mL
Sensitivity: 1 ng/mL
Linear range: 0.1-1000 ng/mL