

**HRP-labeled Anti-p42 ERK2**
RABBIT POLYCLONAL

Catalog Number: AB-231
Quantity: 100 micrograms
Format: PBS
Host: Rabbit
Immunogen: human p42 ERK2 c-terminal peptide-KLH conjugates

Specificity & Preparation: Recognizes p42 MAPK (ERK2), does not cross-react with p44 MAPK (ERK1). Affinity-purified with p42 ERK2 c-terminal peptide conjugated to agarose. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 250 $\mu\text{g/ml}$.

Usage: Reported to be effective for ELISA (0.25 $\mu\text{g/mL}$) and immunoblotting (western blot, 1 $\mu\text{g/mL}$).

Storage: Store the antibody at 4°C. Stable for up to six months. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols