

**Antibody to p42 ERK2 (p42 MAPK)**
RABBIT POLYCLONAL

Catalog Number: AB-232
Quantity: 100 micrograms
Format: TBS, 0.1% NaN₃
Host: Rabbit
Immunogen: human p42 ERK2 c-terminal peptide-KLH conjugates

Specificity & Preparation: Recognizes p42 ERK2 of human rat, bovine, mouse, and rabbit. Does not cross-react with p44 ERK1. Affinity-purified with p42 ERK2 c-terminal peptide conjugated to agarose. The concentration is 250 µg/ml.

Usage: Reported to be effective for ELISA (0.5 µg/ml), immunoblotting (western blot, 1 µg/ml), immunoprecipitation (5 µg/sample).

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols