



6352 CORTE DEL ABETO, STE B  
CARLSBAD, CA 92011 USA  
01.858.642.1988 • WWW.ATSBIO.COM

**Antibody to p70 S6 Kinase**  
Rabbit Polyclonal

**Catalog Number:** AB-241  
**Quantity:** 100 micrograms  
**Format:** TBS, 0.1% NaN<sub>3</sub>  
**Host:** Rabbit  
**Immunogen:** human p70 S6K c-terminal peptide (ISKRPEHLRMNL)-KLH conjugates

**Specificity and Preparation:**

Recognizes p70 S6K of human, mouse, and rat. Affinity-purified with peptide (ISKRPEHLRMNL) conjugated to agarose. The concentration is 250  $\mu\text{g/ml}$ .

**Usage and Storage:**

Reported to be effective for ELISA (0.5  $\mu\text{g/ml}$ ), immunoblotting (western blot, 1  $\mu\text{g/ml}$ ), and immunoprecipitation (5  $\mu\text{g/sample}$ ). Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

**To view protocol(s) for this and other products please visit: [www.ATSBio.com/support/protocols](http://www.ATSBio.com/support/protocols)**