



Antibody to p90 RSK1
RABBIT POLYCLONAL


Catalog Number: AB-244
Quantity: 100 micrograms
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: human p90 RSK1 c-terminal peptide (C-LAQRRVRKLPSTTL)-KLH conjugates

Specificity & Preparation: Recognizes p90 RSK1 of human, mouse, and rat. Affinity-purified with peptide (C-LAQRRVRKLPSTTL) on agarose. The concentration is 250 μ g/ml.

Usage: Reported to be effective for ELISA, immunoblotting (western blot), and immunoprecipitation.

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols



Western blot analysis of p90 RSK1 immunoprecipitated from the mouse brain extract (right) and using RSK1-HRP as a probe. Immunoprecipitation negative control (left).