



6352 CORTE DEL ABETO, STE B  
CARLSBAD, CA 92011 USA  
01.858.642.1988 • WWW.ATSBIO.COM

**HRP-labeled Anti-Phospho-CREBtide**  
Rabbit Polyclonal

**Catalog Number:** AB-246  
**Quantity:** 50 micrograms  
**Format:** PBS, pH 7 with 50% glycerol  
**Host:** Rabbit  
**Immunogen:** phospho-CREBtide (KRREILSRRPpSYR)-KLH

**Specificity and Preparation:**

Recognizes only the phosphorylated CREBtide product of PKA and PKC in a non-radioactive kinase assay, does not cross-react with the non-phosphorylated CREBtide. Affinity-purified with KRREILSRRPpSYR conjugated to agarose followed by affinity adsorption with non-phosphorylated CREBtide. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 250  $\mu$ g/ml.

**Usage and Storage:**

Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:1000). Store the antibody at -20°C. Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

**To view protocol(s) for this and other products please visit: [www.ATSBio.com/support/protocols](http://www.ATSBio.com/support/protocols)**