



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Anti-Phospho-Crosstide
Rabbit Polyclonal

Catalog Number: AB-249
Quantity: 100 micrograms
Format: PBS, pH 7, 0.01% NaN₃
Host: Rabbit
Immunogen: phospho-Crosstide (PRAApTF)-KLH

Specificity and Preparation:

Recognizes only the phosphorylated Crosstide, the product of SGK, PKB in a non-radioactive kinase assay, does not cross-react with the non-phosphorylated Crosstide. Affinity-purified with PRAApTF conjugated to agarose followed by affinity adsorption with non-phosphorylated Crosstide. The concentration is 250 μ g/ml.

Usage and Storage:

Reported to be effective for ELISA, immunoblotting (western blot), and non-radioactive AKT and SGK kinase assay.

Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols