

**Anti-Phosphoserine**
RABBIT POLYCLONAL

Catalog Number: AB-252
Quantity: 100 micrograms
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: phosphoserine-KLH and phosvitin mixture

Specificity & Preparation: Recognizes proteins phosphorylated on serine residues. It does not cross-react with phosphotyrosine. Affinity-purified with phosphoserine conjugated to agarose. The concentration is 250 μ g/ml.

Usage: Reported to be effective for ELISA (1:250), immunoblotting (western blot, 1:250), immunoprecipitation (10 μ g/mg crude protein), and immunohistochemistry (1:50).

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols