

**HRP-labeled Anti-Phosphothreonine**
RABBIT POLYCLONAL

Catalog Number: AB-254
Quantity: 100 micrograms
Format: PBS with 50% glycerol
Host: Rabbit
Immunogen: phosphothreonine-KLH conjugates and phosvitin

Specificity & Preparation: Pan-specific for phosphothreonine, slightly cross-reacts with phosphoserine, but does not cross-react with phosphotyrosine. Affinity-purified with phosphothreonine immobilized conjugated to agarose. Conjugated to horseradish peroxidase (HRP) via reductive amination. The concentration is 1 mg/ml.

Usage: Reported to be effective for ELISA (1:1000) and immunoblotting (western blot, 1:500).

Storage: Store the antibody at -20°C. Stable for up to twelve months. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols