

Antibody to SCD-1
Mouse Monoclonal

Catalog Number: AB-259
Quantity: 100 micrograms
Format: PBS, pH 7.4 with 0.05% sodium azide. Frozen liquid.
Host: Mouse
Isotype: IgG2b
Clone: CD.E10
Immunogen: full-length recombinant protein

Specificity and Preparation:


Reacts with human and mouse proteins. The concentration is 1 mg/ml.

Usage and Storage:

Reported to be effective for immunoblotting (western blot), ELISA, immunoprecipitation, and immunohistochemistry.

Store the antibody at 4°C for three months or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols



Western blot analysis of SCD1 expression in 293 cells.