

**Antibody to Stat5**
RABBIT POLYCLONAL

Catalog Number: AB-262
Quantity: 100 micrograms
Format: PBS, pH 7 with 50% glycerol
Host: Rabbit
Immunogen: human Stat5 (690712 aa) peptide-KLH conjugates

Specificity & Preparation: Recognizes stat5 in human and mouse. Affinity-purified with peptide-affinity chromatography. The concentration is 250 mg/ml.

Usage: Reported to be effective for ELISA (1:2500), immunoblotting (western blot, 1:2500), and immunohistochemistry (1:100).

Storage: Store the antibody at -20°C. Stable for up to three years. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols



Western blot analysis of the extract from mouse spleen with Anti-Stat5.
Right: Spleen from control mouse
Left: Spleen from mouse treated with LPS