

**Antibody to GAT-1, Affinity-Purified**
RABBIT POLYCLONAL

Catalog Number: AB-N37
Quantity: 50 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4)
Host: Rabbit
Immunogen: a peptide from the GAT-1 extracellular domain

Background: GAT-1 is a sodium-coupled neurotransmitter transporter responsible for moving γ -aminobutyric acid (GABA) across cell membranes. GABA is the predominant inhibitory neurotransmitter in the mammalian central nervous system. GAT-1 is widely distributed in both the central and peripheral nervous systems. GAT-1 and GABA are present in numerous neuronal pathways, some of which are implicated in epilepsy, sleep disorders, neuropathic pain, and attention deficit disorders.

Specificity & Preparation: This antibody recognizes the γ -aminobutyric acid (GABA)-1 transporter, GAT-1, in rat. The antibody is an affinity-purified rabbit polyclonal. The peptide used as an antigen has 100% sequence homology between rat, human, mouse, and bovine GAT-1.

Usage: Applications include targeting (targeting agent in GAT1-SAP, Cat. #IT-32).

Storage: Store the antibody at -20°C for one year. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols