



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

**Anti-Conjugated Homocysteic Acid
RABBIT POLYCLONAL**

Catalog Number: AB-T096
Quantity: 50 microliters
Format: Lyophilized and reconstituted with deionized water / 50% glycerol
Host: Rabbit
Isotype: Polyclonal
Immunogen: Synthetic homocysteic acid conjugated to bovine serum albumin

Specificity and Preparation:

Antiserum previously preabsorbed on protein carriers and purified by ammonium sulfate precipitation. This antibody targets conjugated homocysteic acid. **This antibody does not recognize free homocysteic acid.** Using a conjugate Homocysteic Acid-Glutaraldehyde-BSA, antibody specificity was performed with an ELISA test by competition experiments with the following compounds:

COMPOUND	CROSS REACTIVITY §
Homocysteic acid-G-BSA	1
L-Glutamate-G-BSA	1/>1,000
L-Aspartate-G-BSA	1/>1,000

G = Glutaraldehyde, BSA = Bovine Serum Albumin

§ Homocysteic acid-G-BSA concentration/unconjugated or conjugated amino acids concentration at half displacement.

Usage and Storage:

Applications include ELISA (1/1,000-1/5,000). Fixation of tissue for use with these antibodies should be done with glutaraldehyde. The use of paraformaldehyde in conjunction with glutaraldehyde may improve staining results. Store the antibody at 4°C for one month or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

Available Control(s): Homocysteic Acid-BSA

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols