

**Anti-Conjugated Gibberellic Acid**
RAT POLYCLONAL

Catalog Number: AB-T130
Quantity: 50 microliters
Format: Lyophilized and reconstituted with deionized water / 50% glycerol
Host: Rat
Isotype: IgG
Immunogen: Synthetic Gibberellic Acid conjugated to bovine serum albumin

Specificity & Preparation: Antiserum previously preabsorbed on protein carriers and purified by ammonium sulfate precipitation.

This antibody targets conjugated Gibberellic Acid. **This antibody does not recognize free Gibberellic Acid.** Using a conjugate Gibberellic acid-bovine serum albumin (BSA), antibody affinity was performed with an ELISA test by competition experiments.

Usage: Applications include ELISA (1/5,000-1/10,000).

Storage: Store the antibody at 4°C for one month or -20°C in undiluted aliquots for up to one year. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols