

**Antibody to FAK (5B9)**
MOUSE MONOCLONAL

Catalog Number: AB-V64
Quantity: 100 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative.
Host: Mouse
Isotype: IgG1
Clone: 5B9
Immunogen: Affinity-purified tyrosine phosphoproteins from chick embryo fibroblasts expressing activated pp60src

Background: Focal adhesion kinase (FAK) is a non-receptor tyrosine kinase localized at focal adhesions, specialized structures that connect the extracellular matrix to the actin cytoskeleton. Tensin, another focal adhesion component, is a multi-domain protein implicated in cell migration, fibrillar adhesion formation, cartilage development, and the integration of signaling pathways with the cytoskeleton. Tensin has also been identified as a potential therapeutic target in renal disease, wound healing, and cancer.

Specificity & Preparation: This mouse IgG1 antibody is generated against tyrosine phosphoproteins purified from chicken embryonic fibroblasts (CEFs) expressing pp60/c-src 527F and 2A/527F and recognizes chicken FAK.

Usage: Applications include western blot. Working dilutions must be determined by end user.

Storage: Store antibody at -20°C for one year. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

**Selected References:**

1. Kanner SB, Reynolds AB, Vines RR, Parsons JT (1990) Monoclonal antibodies to individual tyrosine-phosphorylated protein substrates of oncogene-encoded tyrosine kinases. Proc Natl Acad Sci U S A 87(9):3328-3332. doi: 10.1073/pnas.87.9.3328 PMID: 2110361

Scan to view
all product
references.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols